

PATCH TESTING : PATIENT TEST RESULTS

Patient Name: _____

ID#: _____

Panel Placement: _____ (date/time)

Panel Removal: _____ (date/time)

Read Date #1: _____

Read Date #2: _____

Read Date #3: _____

Read Date #4: _____

C-1000 Cosmetic Series			Readings				
#	Article#	Description	#1	#2	#3	#4	Comments
1	I-003	Isopropyl myristate 20.0 pet					
2	A-004	Amerchol L 101 50.0 pet					
3	T-016	Triethanolamine 2.0 pet					
4	P-013	Polysorbate 80 / (Polyoxyethylenesorbitanmonooleate(Tween 80)) 5.0 pet					
5	S-004	Sorbitan oleate / (Sorbitan monooleate (Span 80)) 5.0 pet					
6	B-022	2-tert-Butyl-4-methoxyphenol (BHA) 2.0 pet					
7	D-006	BHT / (2,6-Di-tert-butyl-4-cresol) 2.0 pet					
8	O-002	Octyl gallate 0.25 pet					
9	T-014	Triclosan (Irgasan DP300) 2.0 pet					
10	S-003	Sorbic acid 2.0 pet					
11	C-008	p-Chloro-m-cresol / (4-Chloro-3-cresol (PCMC)) 1.0 pet					
12	C-010A	Chloroxylenol (PCMX) / (4-Chloro-3,5-xyleneol (PCMX)) 0.5 pet					
13	T-007	Thimerosal (Merthiolate) 0.1 pet					
14	I-001A	Imidazolidinyl urea (Germall 115) 2.0 pet					
15	H-003	Methenamine / (Hexamethylene tetramine) 2.0 pet					
16	C-005	Chlorhexidine digluconate 0.5 aq					
17	MX-03C	Paraben Mix 16.0 pet					
18	P-008	Phenylmercuric acetate 0.01 aq					
19	C-006	Chloroacetamide / (2-Chloroacetamide) 0.2 pet					
20	H-002	Hexahydro-1,3,5-tris-(2-hydroxyethyl)triazine(Grotan BK) 1.0 aq					
21	C-015	Clioquinol (Chinoform, Vioform) 5.0 pet					
22	E-005	Ethylenediamine dihydrochloride 1.0 pet					
23	A-002	Hydroabietyl alcohol / (Abitol) 10.0 pet					
24	P-011	Phenyl salicylate(Salol) 1.0 pet					
25	H-014C	Benzophenone-3 / (2-Hydroxy-4-methoxybenzophenone) 10.0 pet					
26	S-005	Sorbitan sesquioleate 20.0 pet					
27	P-019A	Propylene glycol 5.0 pet					
28	S-006	Stearyl alcohol 30.0 pet					
29	C-003	Cetyl alcohol 5.0 pet					
30	B-010B	Benzyl salicylate 10.0 pet					
31	B-015A	2-Bromo-2-nitropropane-1,3-diol (Bronopol) 0.25 pet					

PATCH TESTING : PATIENT TEST RESULTS

Patient Name: _____

ID#: _____

Panel Placement: _____ (date/time)

Panel Removal: _____ (date/time)

Read Date #1: _____

Read Date #2: _____

Read Date #3: _____

Read Date #4: _____

C-1000 Cosmetic Series			Readings				Comments
#	Article#	Description	#1	#2	#3	#4	
32	S-002	Sodium-2-pyridinethiol-1-oxide(Sodiumomadine) 0.1 aq					
33	C-018	Cocamidopropylbetaine 1.0 aq					
34	B-008B	Benzyl alcohol 10.0 soft					
35	C-009B	Methylisothiazolinone + Methylchloroisothiazolinone / (Cl+-Me-isothiazolinone (Kathon CG 200ppm)) 0.02 aq					
36	B-028	tert-Butylhydroquinone 1.0 pet					
37	H-016	Drometrizole / (2(2-Hydroxy-5-methyl-phenyl)benzotriazol) 1.0 pet					
38	P-021	Propyl gallate 1.0 pet					
39	D-042	Dodecyl gallate 0.25 pet					
40	C-007A	Quaternium-15 (Dowicil 200) / (1-(3-Chloroallyl)-3,5,7-triaza-1-azoniaadamantane chloride) 1.0 pet					
41	P-025	2-Phenoxyethanol 1.0 pet					
42	D-044A	Diazolidinylurea (Germall II) 2.0 pet					
43	T-036	Tocopherol/(DL-alpha Tocopherol) 100%					
44	D-047A	DMDM Hydantoin 2.0 aq					
45	D-049E	Methyldibromoglutaronitrile (MDBGN) 0.5 pet					
46	T-035B	Tea Tree Oil 5.0 pet					
47	I-008C	Iodopropynyl butyl carbamate 0.2 pet					
48	D-053	Dimethylaminopropylamine 1.0 aq					
49	L-004	Lauryl glucoside 3.0 pet					
50	P-036	Peppermint oil 2.0 pet					
51	S-015	Shellac 20.0 alc					
52	T-037B	Tocopheryl acetate 10.0 pet					
53	T-024B	Turpentine peroxides 0.3 pet					
54	M-035	Methylisothiazolinone 0.02 aq					
55	MX-10B	Musk mix 3.0 pet					
56	O-005	Oleamidopropyl dimethylamine 0.1 aq					