

PATCH TESTING : PATIENT TEST RESULTS

CHEMOTECHNIQUE® DIAGNOSTICS

First in Patch Testing... Since 1981

DORMER LABORATORIES INC.

www.dormer.com

Patient Name: _____

ID#: _____

Panel Placement: _____ (date/time)

Panel Removal: _____ (date/time)

Read Date #1: _____

Read Date #2: _____

Read Date #3: _____

Read Date #4: _____

ME-1000 Medicament Series			Readings				
#	Article#	Description	#1	#2	#3	#4	Comments
1	C-032	Chloramphenicol 5.0 pet					
2	K-001	Kanamycin sulfate 10.0 pet					
3	Q-001	Quinine sulfate 1.0 pet					
4	S-010	Sulfanilamide 5.0 pet					
5	G-006	Gentamicin sulfate 20.0 pet					
6	N-005	Nitrofurazone 1.0 pet					
7	B-032A	Bacitracin 5.0 pet					
8	F-005	Framycetin sulfate 20.0 pet					
9	MX-19	Caine mix III 10.0 pet					
10	M-027	Miconazole 1.0 alc					
11	E-021	Econazole nitrate 1.0 alc					
12	MX-20	Caine mix IV 10.0 pet					
13	F-003	Fusidic acid sodium salt 2.0 pet					
14	T-034	Tioconazole 1.0 pet					