



PATCH TEST RECORD FORM

**MA-1000 (Meth) Acrylate Series Adhesives, Dental & Other
CHEMOTECHNIQUE PATCH TEST ALLERGENS**

Patient Name: _____

Date Applied: _____

Date of 1st reading: _____

Date of 2nd reading: _____

Date of 3rd reading: _____

Compound	Conc %	48H	96H	7 days
1 Methyl methacrylate	2.0 pet	[]	[]	[]
2 n-Butyl methacrylate	2.0 pet	[]	[]	[]
3 2-Hydroxyethyl methacrylate	2.0 pet	[]	[]	[]
4 2-Hydroxypropyl methacrylate	2.0 pet	[]	[]	[]
5 Ethyleneglycol dimethacrylate	2.0 pet	[]	[]	[]
6 Triethyleneglycol dimethacrylate	2.0 pet	[]	[]	[]
7 1,4-Butanediol dimethacrylate	2.0 pet	[]	[]	[]
8 Urethane dimethacrylate	2.0 pet	[]	[]	[]
9 BIS-MA	2.0 pet	[]	[]	[]
10 BIS-GMA	2.0 pet	[]	[]	[]
11 1,6-Hexandiol diacrylate	0.1 pet	[]	[]	[]
12 Tetrahydrofurfuryl methacrylate	2.0 pet	[]	[]	[]
13 Tetraethyleneglycol dimethacrylate	2.0 pet	[]	[]	[]
14 N,N-Dimethylaminoethyl methacrylate	0.2 pet	[]	[]	[]
15 Ethyl cyanoacrylate	10.0 pet	[]	[]	[]