



PATCH TEST RECORD FORM
P-1000 Photographic Chemicals Series
CHEMOTECHNIQUE PATCH TEST ALLERGENS

Patient Name: _____

Date Applied: _____

Date of 1st reading: _____

Date of 2nd reading: _____

Date of 3rd reading: _____

	Compound	Conc %	48H	96H	7 days
1	CD-2, N,N-Diethyl-2-methyl-1,4-phenylenediamine-HCl	1.0 pet	[]	[]	[]
2	CD-3, n-Ethyl-N-(2-methanesulfonamidoethyl)-2-methyl	1.0 pet	[]	[]	[]
3	CD-4, N-ethyl-N-(2-hydroxy-ethyl)-2-methyl-1,4-PPD-sulfate	1.0 pet	[]	[]	[]
4	4-Methylaminophenol sulfate (Metol)	1.0 pet	[]	[]	[]
5	Hydroquinone	1.0 pet	[]	[]	[]
6	Phenidone	1.0 pet	[]	[]	[]
7	Hydroxylammonium chloride	0.1 aq	[]	[]	[]
8	Ammonium persulfate	2.5 pet	[]	[]	[]
9	Ethylenediamine dihydrochloride	1.0 pet	[]	[]	[]
10	1H-Benzotriazole	1.0 pet	[]	[]	[]
11	Glutaraldehyde	0.2 pet	[]	[]	[]
12	Benzyl alcohol	10.0 soft	[]	[]	[]
13	Hydroxylammonium sulfate	0.1 aq	[]	[]	[]
14	Potassium dichromate	0.5 pet	[]	[]	[]
15	4-Amino-N,N-diethylaniline sulfate (TSS)	1.0 pet	[]	[]	[]
16	Tricresyl phosphate	5.0 pet	[]	[]	[]